

## SHAPE CAPTURE MARK-UP ● ● ● ●

Use these standardized mark-ups directly on the shape capture bag to improve accuracy and consistency. Clear markings streamline the transition from capture to molding, reduce field adjustments, and support faster delivery of a precise final product.

**HR**

### Headrest placement\*

*\*Marked on front; mounting bracket placement will be transferred to rear surface*

**+ # " ↑**

### Add back support height

*+2" ↑ Add 2 inches of back height*

**(R)**

### Add Digital Relief or Digital Offload area

*At a Gibbus, Rib Hump, or for cushion offloading*

**(X)**

### Add T-nuts

*For Chest Strap or Harness*

**+ # " ← or + # " →**

### Increase width of lateral

*+1.5" → Extend lateral by 1.5 inches*

**N**

### Add Notches (seat belt, back canes, etc)

*Indicate appropriate location and shape*

**PRZ**

### Add Pressure Relief Zone (soft spot)

*Outline with appropriate shape*

**+ # " ↓**

### Add depth for leg length discrepancy

*+1" ↓ Add 1 inch of depth at leg trough*

**+ # " ⤿**

### Abductor/Pommel location and size

*+1" Add 1 inch of height/cushion thickness*

**+ # " ↑**

### Add height at lateral thigh area

*+2" ↑ Add 2 inches of lateral thigh height*