

Order Form

www.matrixseatingUSA.com info@matrixseatingusa.com (800) 986-9319

	ITEN	AIZED	
Please note: For Original Matrix order	forms, reach us a	t (800) 986-9319 or email info@matrixseatingusa.	com
BACKS		CUSHIONS	
AIREZ Custom Back	(MACB)	AIREZ Custom Seat Cushion	(MACSC)
AIREZ Custom Back + ABS-X Shell	(MACB) + (MAB-ABS)	Auxafoam Custom Seat Cushion	(MCSC)
PROVIDER INFORMATIO Account #:			
Company Name:		Contact Name:	
Phone Number:		Email:	
Mailing Address:			
City:		State: ZIP:	
Shipping Address:			
City:		State: ZIP:	
Our seating products must be installed by Ma	trix Certified Provid	lers and will not be drop shipped to end users.	
CLIENT INFORMATION			
Client First and Last Name:			
Client Height:		Client Weight*:	
Wheelchair Make:		Model:	
Frame Width:		Frame Depth:	
Seat Base Type: O Solid Pan	Sling	Seat	
Special Considerations*: D Hi Bold sections indicate required fields. * Upgrades may be required for clients whose	-	 Spasticity 50 LBS or present with severe spasticity or extensor tone. 	\bigcirc
MEASUREMENTS		ſ	25
Shoulder Width:	Buttocks to	Top of Shoulders**:	1
Hip Width:	Seat Base to	b Back of Knee, L:	
** Additional hardware may be required if mo		Back of Knee, R:	
BE	SURE TO A	LSO INCLUDE	

• Photos of the front and side of the client during shape capture.

• Photos of the captured shape.

Custom Backs

• AIREZ

Matrix AIREZ Custom Back (MACB)

\$2175

For maximum airflow and the lightest weight with strength and rigidity.

Standard configuration includes: an ultra-lightweight adjustable back, lightweight + breathable 3D mesh liner, and a breathable spacer fabric cover.

SPECIFY BACK WIDTH (SELECT ONE)

- Standard Trochanter Width <20"</p>
- □ XL Trochanter Width 21"-24"

HARDWARE OPTIONS (SELECT ONE)

Ultra lightweight mounting hardware, mounting width:

– 🗆 Frame width: _____"

Q Radian 360 mounting hardware, mounting width:

- Frame width: _____"
- D Add additional Radian 360 of same size*
- Fixed Pivot mounting hardware, system size**:
 - D Small Medium, 4 mounting points
 - D Large, 6 mounting points
 - 🗆 Extra Large, 8 mounting points
- Power Wheelchair Adapter***
 For mounting on a solid back.
- Power Wheelchair Bundle*** (PWB)
 For use with wheelchair manufacture's back interface. Includes surface mount hardware, lightweight hardware, and straight bar.

Omit Hardware
(RAD-O)



Power Wheelchair Adapter (PWA)



Power Wheelchair Hardware Bundle (PWB)

* Required for clients whose back is equal to or taller than 29" or whose weight exceeds 250 lbs or who presents with severe spasticity or extensor tone. ** Requires ABS-X outer shell.

*** Requires 3" of seat depth adjustment to accommodate hardware.

Page 2/9



OAIREZ + ABS-X

Matrix AIREZ Custom

Back (MACB) + (MAB-ABS)

* Required for clients whose back is equal to or taller than 29" or

whose weight exceeds 250 lbs or

\$395

\$595

\$595

\$866.40

\$1299.60

\$1732.80

\$645

\$495

\$2175 + \$485 The same AIREZ back with an

additional ABS-X shell*.

who presents with severe

spasticity or extensor tone.

(MAB-TSW)

(VTMH)

(RAD)

(RAD-EX)

(PIV-4)

(PIV-6)

(PIV-8)

(PWA)

(MAB-TXLW) \$347

Ultra Lightweight Mounting Hardware (VTMH)



Radian 360 Hardware (RAD)



Fixed Pivot Hardware (PIV)

Custom Backs (E2617) + Custom Seat Cushions (E2609

Published Jun 8, 2023

info@matrixseatingusa.com

(800) 986-9319



ITEMIZED

(E2617)

HARDWARE MOUNTING ACCESSORIES

Cane Claw Hardware - Round Mount	(R360-CC) No Charge	
O Surface Mount Hardware - Flat Mount	(R360-SM) No Charge	Cane Claw Hardware
Quick Release Hardware Upgrade	(R360-QR) \$95	(R360-CC)
Additional Pair of Mounting Clamps		
🖵 🗆 Cane Claw	(R360-EC) \$230	
🖵 🛛 Surface	(R360-ES) \$230	Surface Mount for Non-cane Attachment (R360-SM)

(AFBU)

\$177

OPTIONAL - LINER UPGRADE

O Therapeutic Premium 3D Liner*

* Recommended for users with compromised skin integrity and/or requiring tilt for pressure relief.

OPTIONAL - ALTERNATIVE BACK COVER

Replace with Dartex Style Cover	(CVRBU-DAR)	\$205
---------------------------------	-------------	-------

Dartex Style Covers are silver infused for anti-microbial properties and smoother skin contact.

OPTIONAL - ADDITIONAL BACK COVERS

(CVRB-SF)	\$345
(CVRB-DAR)	\$550
(CVRB-IP)	\$250
(MCB-PF)	\$150
	(CVRB-DAR) (CVRB-IP)



Spacer Fabric Cover

* Only recommended for patients with chronic incontinence.

OPTIONAL - BACK ACCESSORIES

O Universal Headrest Mount	
– 🗆 Omit installation	
C p Accessory Mounting Rails	
🖵 🛛 Shoulder Harness Guides	
🖵 🗆 Omit installation	

н

UHM with Accessory Mounting Rails

(CFI-HR) + (MCB-IM)

(CFI-HR)	\$175	
(CFI-O)	No Charge	
(MCB-IM)	\$271	
(MCB-SHG)	\$200	
(MCB-O)	No Charge	



Universal Headrest Mount (CFI-HR)



UHM with Accessory Mounting Rails and Shoulder Harness Guides (CFI-HR) + (MCB-IM) + (MCB-SHG)



info@matrixseatingusa.com

Published Jun 8, 2023

ITEMIZED

Custom Backs (E2617) + Custom Seat Cushions (E2609)

OPTIONAL - BACK SPECIFICATIONS

If no specifications are selected, the following will be manufactured as captured.

Ο	Extend patient LEFT lateral thoracic support" forward of line	(MCB-EDL-L)	\$290
Ο	Extend patient RIGHT lateral thoracic support" forward of line	(MCB-EDL-R)	\$290
Ο	Increase patient LEFT lateral thoracic support" above line	(MCB-IHL-L)	\$190
Ο	Increase patient RIGHT lateral thoracic support" above line	(MCB-IHL-R)	\$190
Ο	Extend back height" above line	(MCB-EBH)	\$290

Cushion products available on following page

ITEMIZED

info@matrixseatingusa.com

Published Jun 8, 2023

Custom Cushions		JXAF		2609)	Cus
					ustom Backs
Matrix AIREZ Custom Seat Cushior \$2065 Includes a custom-molded AIREZ cushion and brea fabric cover.	athable spacer Includes	a custom-molded Au	t Cushion (MCSC) \$2065 xafoam cushion, breathable space l incontinence-proof interior cover.		i (E2617) + C
WHEELCHAIR SEAT SPECI Seat Base Type: O Solid Pan	-	T ONE)			+ Cu
SPECIFY SEAT WIDTH (SEI	LECT ONE) (SS-1	AIRI N)	EZ AUXAFOA	АМ	ustom S
 XL Width 21" - 24" Specify Exact Width:" Tapered Width Front Width:" Back Width 	(SS-XL (SC-T				eat C
 OPTIONAL - ALTERNATIVE Replace with Dartex Style Cover Dartex Style Covers are silver infused for and 	er (CVRSU				ushions (E2
OPTIONAL - ADDITIONAL O Spacer Fabric	SEAT COVERS	- SF) \$200			2609)
 Dartex Style (Smoother Skin Cor Incontinence-proof interior cov Incontinence-proof exterior cov 	ver with zip* (CVRS	I-IP) \$250		/	
* Not recommended for Matrix AIREZ Cushic ** Only recommended for patients with chro	ons.	,			
 OPTIONAL - SEAT ACCESS Integrated Sling Seat Mounting *** Recommended for chairs with sling seats 	g Platform (SSM				ITEMIZED
Addition	nal seat options available on t	ollowing page			E D
Page 5/9 www.matrixseatingUSA.com	info@matrixseatir		Published	Jun 8, 2023 986-9319	Æ

OPTIONAL - SEAT SPECIFICATIONS

If no specifications are selected, the following will be manufactured as captured.

Modify Seating Height Increase overall height by" As low as possible 	(SC-THI) (SC-THD)	\$140 \$140		tom Backs
Modify Length - PATIENT LEFT Longer than shape capture by" Shorter than shape capture by" 	(SC-CL-ALL) (SC-CL-ALS)	\$140 \$140		ks (E2617)
Modify Length - PATIENT RIGHT Longer than shape capture by" Shorter than shape capture by" 	(SC-CL-ARL) (SC-CL-ARS)	\$140 \$140		+ Custom
 Modify Medial Thigh Support □ Eliminate completely □ Increase height by" □ Decrease height by" 	(SC-MTE) (SC-MTI) (SC-MTD)	No Charge \$110 \$110		Seat Cushions
 Modify Lateral Thigh Support - PATIENT LEFT Eliminate completely Increase height by" Decrease height by" 	(SC-LLTE) (SC-LLTI) (SC-LLTD)	No Charge \$110 No Charge		ons (E2609)
Modify Lateral Thigh Support - PATIENT RIGHT Eliminate completely Increase height by" Decrease height by" 		No Charge \$110 No Charge		
 Lateral Thigh Support Reinforcement* * Increases lateral thigh support rigidity. Applicable to Auxafoam of Soft Spot Adjustment Must notate location on shape capture. Seat Positioning Tray Fitted tray that allows for repeatable placement of cushion to ensign 	(SC-SS) (SC-SPT)	\$190 \$125 \$350 hing.	Seat Positioning Tray (SC-SPT)	ITEMIZED

Page 6/9

info@matrixseatingusa.com

Published Jun 8, 2023 (800) 986-9319

Thank you [•]	for choosing	Matrix	Seating	USA
------------------------	--------------	--------	---------	-----

ADDITIONAL COMMENTS OR INSTRUCTIONS

ITEMIZED

 \mathcal{C}

Page 7/9

Published Jun 8, 2023

Hardware Guide

ULTRA LIGHTWEIGHT MOUNTING HARDWARE

 The Matrix Ultra lightweight hardware is the lightest hardware Matrix Seating USA offers. While extremely lightweight and low profile, it has limited adjustability to and from the client and limited rotation on the X and Y axes. Height adjustment is achieved though adjusting the clamp receivers on the canes. It is typically used for mounting to self-propelled manual wheelchairs.

RADIAN 360

• Matrix Radian 360 offers the greatest level of adjustability and strength. As the name suggests, the Radian 360's ball and socket configuration enables an infinite level of adjustability. It is typically used for mounting to tilt in space wheelchairs and power wheelchairs with available canes.

FIXED PIVOT HARDWARE

• Fixed pivot hardware is available for mounting nontraditional seating components. Examples include: seating systems that have the back and seat combined, supine positioning systems mounted to a tilt in space wheelchair, or creative mounting solutions for backs with extreme curvature. Fixed pivot hardware requires the ABS-X outer shell upgrade.

POWER WHEELCHAIR ADAPTER

• The Power Wheelchair Adapter (PWA) is the primary solution for mounting to powered wheelchairs with solid backs 7" or wider via (4) bolts. The PWA allows for rotation and X/Y axis adjustability to enable easy client positioning after initial installation while using under 3" of seat depth.

POWER WHEELCHAIR BUNDLE

 The Power Wheelchair Bundle (PWB) is utilized for mounting to powered wheelchairs without canes or with solid backs with cut-outs that does not enable the Powered Wheelchair Adapter to be utilized. Paired with the VTMH-S straight horizontal mounting shaft and R360-SM surface mount receivers, the hardware is mounted on available surfaces 10" or greater between the canes via (8) bolts, (4) in each surface mount receiver using under 2.5" of seat depth.





ITEMIZED





www.matrixseatingUSA.com

info@matrixseatingusa.com



Recommendations by wheelchair type

MANUAL WHEELCHAIRS

AIREZ Back

- Lightest weight option with full airflow benefits.
- Ultra Lightweight Hardware
 - Lightest mounting solution with appropriate adjustability for most use cases.
- AIREZ Cushion
 - Lightest weight option with superior airflow benefits and increased user comfort.
- Spacer Fabric Covers
 - Enables the greatest airflow benefits AIREZ has to offer.

POWER WHEELCHAIRS (PWC)

- Options based on wheelchair mounting type:
- -• PWC w/ canes; Radian 360
 - No seat depth adjustment required.
- • PWC w/ solid back; Power Wheelchair Adapter
 - Increases require 3" seat depth to accommodate hardware.
- Cover selection depending on cushion type:
- • AIREZ; Spacer Fabric Cover
 - Enables the greatest airflow benefits AIREZ has to offer.
- Auxafoam Seat; Dartex or Spacer Fabric with fluid proof liner
 - Fluid proof liner protect foam cushion from moisture.

TILT IN SPACE WHEELCHAIRS

- AIREZ Back + Therapeutic Premium 3D liner
- ABS-X outer shell
 - Increases levels of inner shell shape adjustability, rigidity for extreme tone/spasticity and adds easily accessible mounting
 points for chest straps and other accessories.
- Radian 360 Hardware
 - Greater adjustability options for wheelchair mounting.
- Does the client possess sensitive skin?
 - Yes; Dartex Style cover (smooth)
 - Especially useful for clients who are diaper-only with leg exposure to seating surface.
 - No; Spacer Fabric cover
 - Helps to maintain patient positioning and AIREZ's airflow benefits.

ITEMIZED

Published Jun 8, 2023

Page 9/9